


## THE WORD ACCURACY GETS USED A LOT.... BUT WHAT DOES ACCURACY REALLY MEAN?

“In science and engineering, accuracy is the degree of conformity of a measured or calculated value to its actual or specified value. Compare and contrast with precision” -Websters Dictionary.

Based on this definition, all of the following meters are accurate, because they all fall within the manufacturer’s accepted range. The question is:

Is the manufacturer’s acceptable range good enough for your patients?

### NORMAL CONTROL RANGE\*



NORMAL CONTROL RANGE		NORMAL CONTROL RANGE	
BIONIME GM100	82 - 110	HDI TRUE TRACK	176 - 238
ABBOTT FREESTYLE LITE	83 - 125	WAVESENSE PRESTO	110 - 157
BAYER CONTOUR	106 - 147	ARKRAY GLUCOCARD	97 - 162
ROCHE	42 - 72		

\*Based on products purchased at retail May 2009.

<b>Job Title:</b> Bionime Control Solution Range Sell Sheet   <b>Job #</b> 188-1101			<b>SAGON   PHIOR:</b> 310.575.4441		<b>MECH</b>
<b>Date:</b> 5.28.09	<b>Media/Pub:</b>	<b>Size: Trim:</b> 8.5" x 11"			
<b>Project Mgr:</b> LJ Feinstein	<b>Creative:</b> Bill Fate	<b>Notes: Color:</b> Four Color Process			

### Strip Technology

Biomime strips are proven reliable and accurate. Our patented, noble-metal nano-technology delivers laboratory quality results. They're also very easy to use. We made them larger, gave them a textured grip, and craft them out of rigid plastic instead of the flimsy paper or Mylar other manufacturers use. They're perfect for people who have lost mobility or vision as a result of complications from advanced diabetes, as well as for anyone who appreciates quality design and construction.



**GS 300**

Noble Gold Metal Electrode Test Strip

- Patented design eliminates contact with sample
- Noble gold technology for increased accuracy
- Ergonomic and easy to handle



**GS 100**

### Lancet Device

The GD500 is a reusable lancing device for home and professional use. With state-of-the-art technology, the GD500 automatically injects and retracts a lancet, reducing pain and preventing contamination. Safe and flexible, the GD500 allows for testing from the finger, palm or forearm, and its seven-depth adjustable cap allows individual users to choose the setting most comfortable and effective for them. Portable and practical, the GD500 makes testing convenient, anywhere, anytime.



**GD 500**

- Award-winning design
- Special safety switch
- Lancet eject function
- Engineered to reduce pain
- 7 depth selections
- Compatible with most lancets
- Alternate site testing (finger, palm or forearm)

## BIONIME

Accurate intel for the war on diabetes.

#### BIONIME USA CORPORATION

11760 Sorrento Valley Road  
Suite G San Diego  
CA 92121, USA  
Tel: +1 858 481-8485  
Fax: +1 858 481-8472  
www.bionimeusa.com

#### BIONIME CORPORATION

694 Renhua Road, Dali City  
Taichung County, Taiwan 412  
www.bionime.com

#### BIONIME GMBH

Tramstrasse 16  
9442  
Berneck, Switzerland  
www.bionime.ch

#### BIONIME (SHENZHEN) CO., LTD.

No 708, 7F Beauty AAA Building  
North Ren Min Rd, Longhua  
Bao'an District, ShenZhen City, China

Copyright © 2009 BIONIME www.bionimeusa.com

<b>Job Title:</b> Bionime Control Solution Range Sell Sheet   <b>Job #</b> 188-1101		<b>SAGON   PHIOR:</b> 310.575.4441		<b>MECH</b>
<b>Date:</b> 5.28.09	<b>Media/Pub:</b>	<b>Size: Trim:</b> 8.5" x 11"		
<b>Project Mgr:</b> LJ Feinstein	<b>Creative:</b> Bill Fate	<b>Notes: Color:</b> Four Color Process		